



# Light

## Next Generation Science Standards

### NGSS Science and Engineering Practices:

- Asking questions and defining problems
- Developing and using models
- Planning and carrying out investigations
- Analyzing and interpreting data
- Using mathematics and computational thinking
- Constructing explanations and designing solutions
- Engaging in argument from evidence
- Obtaining, evaluating, and communicating information

### NGSS Cross-cutting Concepts:

- Patterns
- Cause and effect
- Scale, proportion, and quantity
- Systems and system models
- Energy and matter
- Structure and function
- Stability and change

### NGSS Disciplinary Core Ideas:

- PS3.D: Energy in Chemical Processes and Everyday Life

## Initial Prep Time

Approx. 10-15 min. per apparatus

## Lesson Time

1 – 2 class periods, depending on experiments completed

## Assembly Requirements

- Scissors
- Small Philips screwdriver

### Materials (for each lab group):

- Horizon Solar Hydrogen Science Kit
- Distilled water
- AA batteries
- Protractor
- Stopwatch
- Colored construction paper
- Various colored light filters
- Heat lamp and/or UV lamp (optional)
- Horizon Renewable Energy Monitor or multimeter (optional)



# Light



## Lab Setup

- We recommend completing steps 1 and 2 in the Assembly Guide for each electrolyzer so your students do not have to cut tubing or fill the electrolyzer initially.
- Lab includes small parts that can go missing easily. Set up a resource area for each lab table or for the entire class to minimize lost pieces.
- A heat lamp or UV lamp may be used during experiment #2, if available.
- If you don't have access to a multimeter or Horizon Renewable Energy Monitor, omit the Measurements section of this activity.



## Safety

- Battery packs can short out and heat up if the red and black contacts touch each other while the unit is in the on position. Be sure to keep them off when not in use.
- Using regular tap water instead of distilled water will severely shorten the lifespan of the fuel cells. Distilled water can be found at most pharmacies or drug stores.
- Running electric current through dry fuel cells or attaching the battery packs backwards can destroy the fuel cells. Be sure to always connect red to red and black to black.
- Beware of water spills, and don't be surprised if someone tries to start a syringe water fight.



## Notes on the Solar Hydrogen Kit

- Direct sunlight, or a strong electric light, is necessary for operation. Overcast and indirect sunlight may not provide sufficient energy. Be sure any artificial light source is close to the solar panel.
- Be sure to line up the gaps on the inner cylinders of the H<sub>2</sub> and O<sub>2</sub> tanks so that bubbles can escape.
- You may need to inject more water into the O<sub>2</sub> side of the cell if the electrolysis reaction is being sluggish. Wait 3 minutes and then try again.



## Common Problems

- The motor's fan sometimes needs a little push to get started.
- If there's hydrogen left but the motor doesn't run, you may have to purge the fuel cell. Unplug the black plug and then quickly replace it to purge impure gases.
- If the water level doesn't change after purging the cells, make sure the gaps on the base of the inner cylinders are open so that water can fill them.



# Light



## Goals

- ✓ Use a solar panel to generate electricity from light
- ✓ Run a motor with the electricity generated
- ✓ Use the speed of the motor to measure light energy



## Background

Light is a strange phenomenon. You've probably been using two highly sensitive light detectors since the day you were born, and they're helping you to read these words right now. But what we see as light is just part of a diverse type of energy that exists all over the universe and has many uses here on our own planet as well.

Light is just a small part of something known as the electromagnetic spectrum, a form of energy that travels through space as waves. You can see only part of that spectrum with your eyes, which your brain interprets as colors. Difference in wavelength (the distance between the peaks of the waves) result in different colors. The colors you can see range from red at the long end of the spectrum to violet at the short end.

But there are many more "colors" beyond those that you can't see, although you may have heard of their names. We call the colors with wavelengths too short to see "ultraviolet" and those with wavelengths too long to see "infrared." Other types of electromagnetic waves, like X-rays and gamma rays, have even shorter wavelengths than ultraviolet. Radio waves and microwaves have even longer wavelengths than infrared.

Solar power is a way of generating electricity that uses the energy contained in these waves to create an electric current. During this activity, you'll use a solar panel to generate an electric current and describe how it works.



## Procedure

1. Use your solar cell to power the small motor that controls the fan. You'll need to connect the solar cell to the fan using wires to carry the electricity. Why do you think you need two wires?
2. When you've connected the solar cell to the motor, you may have to give the fan a little push to get it started. The solar cell will work best in direct sunlight. What happens to the fan if you try the solar cell with other light sources?
3. You can use the electricity from the solar panel to generate hydrogen gas using the electrolyzer. The electrolyzer is the square with "H<sub>2</sub>" and "O<sub>2</sub>" printed on either side. What do you think will happen if you connect it to a source of electricity like the solar cell?
4. Your electrolyzer is also a hydrogen fuel cell that can generate electricity from hydrogen and oxygen. It has two small tubes attached to it. Is there anywhere else on the fuel cell that you could attach the longer tubes?
5. Look at the remaining pieces of your kit. If the fuel cell splits water into hydrogen and oxygen gases, what could you use to trap the gases so they don't float away?



# Light

6. Connect the tubes of your fuel cell so that you can trap the gases. To generate hydrogen, you'll need to supply an electric current. You can do this with the battery pack or the solar cell. Try both. Which is better at producing hydrogen? How do you know?
7. When you've produced hydrogen, you can use the fuel cell to power the motor just like you did with the solar cell. Plug the motor into the fuel cell and it should start turning. Note in your observations if you see any difference in how the motor works with the fuel cell instead of the solar cell.
8. When you've produced hydrogen, you can use the fuel cell to power the motor just like you did with the solar cell. Plug the motor into the fuel cell and it should start turning. Note in your observations if you see any difference in how the motor works with the fuel cell instead of the solar cell.



## Observations



## Experimentation

1. You can use colored plastic gels, or different lightbulbs, to change the color of light hitting the solar panel. Do certain colors work better than others? Try using the solar panel to run the electrolyzer while the panel is hit with different wavelengths of light and record your observations below:

Light Color:	Time to fill H <sub>2</sub> :	Observations:



# Light

2. The solar panel is dark in color. Does the color of its surroundings affect how well it collects light for generating electricity? Try using the panel to run the electrolyzer while the panel is in front of different colored backgrounds and record your observations below:

Background Color:	Time to fill H2:	Observations:

3. Attach the solar panel to the motor and use a piece of paper or other method to shade parts of the panel and observe the effects. How much of the solar panel can you cover before the motor stops running? Does it matter which side of the solar panel is shaded?

**Students should note that, depending on which side you shade, it doesn't take much at all to stop the motor. This is the result of how the individual photovoltaic cells in the solar cell are wired together.**



## Measurement

For this section, you will need a multimeter or the Horizon Renewable Energy Monitor. For an introduction to using a multimeter, [click here](#).

1. Attach the motor to your solar panel. Measure the current in Amps and the voltage in Volts while shading the solar panel to find the minimum values for each that will still run the motor. Record your answers below:

**(Answers will vary, but check that they are within reason, i.e. not >1A.)**



# Light

Current: \_\_\_\_\_ A

Voltage: \_\_\_\_\_ V

2. Voltage is equal to the current multiplied by the resistance ( $V = IR$ ), so according to your data what is the resistance of the motor?

Resistance: \_\_\_\_\_  $\Omega$

3. Use different colors of light with your solar panel as before. Measure the current in Amps and the voltage in Volts while running the motor. What color gave the highest values? Record your answers below:

Color: \_\_\_\_\_

Current: \_\_\_\_\_ A

Voltage: \_\_\_\_\_ V



## Analysis

1. Make a scientific claim about what you observed while running the solar panel.

**Claim should reference the limits of the solar cell's capabilities, in terms of wavelengths of light, amount of light, or absorption of its surrounding.**

*Example: "The solar panel works best with visible wavelengths of light."*



## Light

2. What evidence do you have to back up your scientific claim?

**Evidence should cite data in Observations and/or Experimentation sections.**

*Example: "The electrolyzer filled the hydrogen cylinder in 102 seconds when the solar panel was under visible light, but it took 240 seconds under ultraviolet and 307 seconds under infrared."*

3. What reasoning did you use to support your claim?

**Reasoning can draw from Background section and/or other materials used in class.**

*Example: "The solar panel must have been producing less current if it took longer to power the electrolyzer."*

4. Design an experiment that could test the relationship between the energy of light and its wavelength.

**There are many possible answers, but there should be a mention of a way to measure both the wavelength and energy of the light, and clear control and experimental groups in the experiment.**



## Conclusions

1. Based on your observations, do you think a solar panel would be useful for generating electric energy from any type of light? Explain your reasoning.

**"Yes" or "no" are both acceptable answers, so long as students are able to point to specific data from their experiments to back up their assertion.**

2. What would you say is the most important factor in determining how much electric energy a solar panel produces?

**Student answers should reference data collected in all experiments.**



## Light

3. Based on your observations, what color of light emits the most energy?

**Answers will depend on the variety of colors used.**

4. Based on your observations, what color of background absorbs the most energy?

**Answers will depend on the variety of colors used.**